

## GSPT2 抗原（重组蛋白）

中文名称：GSPT2 抗原（重组蛋白）

英文名称：GSPT2 Antigen (Recombinant Protein)

别名：GST2; ERF3B

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 429-628 amino acids of human GSPT2

技术规格：

<b>Full name:</b>	G1 to S phase transition 2
<b>Synonyms:</b>	GST2; ERF3B
<b>Swissprot:</b>	Q8IYD1
<b>Gene Accession:</b>	BC036077
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a GTPase that belongs to the GTP-binding elongation factor family. The encoded protein is a polypeptide release factor that complexes with eukaryotic peptide chain release factor 1 to mediate translation termination. This protein may also be involved in mRNA stability.