

兔抗 FKBP6 多克隆抗体

中文名称：兔抗 FKBP6 多克隆抗体

英文名称：Anti-FKBP6 rabbit polyclonal antibody

别名：FKBP36

相关类别：一抗

储存：冷冻（-20℃）

抗原：FKBP6

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a cis-trans peptidyl-prolyl isomerase that may function in immunoregulation and basic cellular processes involving protein folding and trafficking. This gene is located in a chromosomal region that is deleted in Williams-Beuren syndrome. Defects in this gene may cause male infertility. There are multiple pseudogenes for this gene located nearby on chromosome 7. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	FKBP6

Immunogen:	Full length fusion protein
Full name:	FK506 binding protein 6, 36kDa
Synonyms :	FKBP36
SwissProt:	O75344
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human esophagus cancer
IHC Recommend dilution:	25-100



