

兔抗 FKRP 多克隆抗体

中文名称：兔抗 FKRP 多克隆抗体

英文名称： Anti-FKRP rabbit polyclonal antibody

别名： MDC1C; LGMD2I; MDDGA5; MDDGB5; MDDGC5

抗原： FKRP

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a protein which is targeted to the medial Golgi apparatus and is necessary for posttranslational modification of dystroglycan. Mutations in this gene have been associated with congenital muscular dystrophy, mental retardation, and cerebellar cysts. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined.
Applications:	WB
Name of antibody:	FKRP
Immunogen:	Synthesized peptide derived from N-terminal of huma

	n FKRP.
Full name:	fukutin related protein
Synonyms :	MDC1C; LGMD2I; MDDGA5; MDDGB5; MDDGC5
SwissProt:	Q9H9S5
WB Predicted band size:	55 kDa
WB Positive control:	Jurkat cells and COLO205 lysates
WB Recommended dilution:	500-3000

