

兔抗 FLCN 多克隆抗体

中文名称：兔抗 FLCN 多克隆抗体

英文名称： Anti-FLCN rabbit polyclonal antibody

别名： BHD; FLCL

相关类别： 一抗

抗原： FLCN

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is located within the Smith-Magenis syndrome region on chromosome 17. Mutations in this gene are associated with Birt-Hogg-Dube syndrome, which is characterized by fibrofolliculomas, renal tumors, lung cysts, and pneumothorax. Alternative splicing of this gene results in two transcript variants encoding different isoforms.
Applications:	ELISA, IHC
Name of antibody:	FLCN
Immunogen:	Fusion protein of human FLCN
Full name:	folliculin

Synonyms :	BHD; FLCL
SwissProt:	Q8NFG4
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human thyroid cancer and Human liver cancer
IHC Recommend dilution:	100-300

