

HMGCLL1 抗原（重组蛋白）

中文名称：HMGCLL1 抗原（重组蛋白）

英文名称：HMGCLL1 Antigen (Recombinant Protein)

别名：er-cHL; bA418P12.1

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 171-370 amino acids of human HMGCLL1

技术规格：

Full name:	3-hydroxymethyl-3-methylglutaryl-CoA lyase like 1
Synonyms:	er-cHL; bA418P12.1
Swissprot:	Q8TB92
Gene Accession:	BC024194
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Non-mitochondrial 3-hydroxymethyl-3-methylglutaryl-CoA lyase that catalyzes a cation-dependent cleavage of (S)-3-hydroxy-3-methylglutaryl-CoA into acetyl-CoA and acetoacetate, a key step in ketogenesis, the products of which support energy production in nonhepatic animal tissues.