

HS1BP3 抗原（重组蛋白）

中文名称：HS1BP3 抗原（重组蛋白）

英文名称：HS1BP3 Antigen (Recombinant Protein)

别名：ETM2; HS1-BP3

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 193-392 amino acids of human HS1BP3

技术规格：

Full name:	HCLS1 binding protein 3
Synonyms:	ETM2; HS1-BP3
Swissprot:	Q53T59
Gene Accession:	BC050636
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene shares similarity with mouse Hs1bp3, an Hcls1/Hs1-interacting protein that may be involved in lymphocyte activation.