

## IFT22 抗原（重组蛋白）

中文名称：IFT22 抗原（重组蛋白）

英文名称：IFT22 Antigen (Recombinant Protein)

别名：RABL5

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 40-185 amino acids of human IFT22

技术规格：

<b>Full name:</b>	intraflagellar transport 22
<b>Synonyms:</b>	RABL5
<b>Swissprot:</b>	Q9H7X7
<b>Gene Accession:</b>	BC009823
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Small GTPase-like component of the intraflagellar transport (IFT) complex B.