

## IFT27 抗原（重组蛋白）

中文名称：IFT27 抗原（重组蛋白）

英文名称：IFT27 Antigen (Recombinant Protein)

别名：RAYL; BBS19; RABL4

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-186 amino acids of human IFT27

技术规格：

<b>Full name:</b>	intraflagellar transport 27
<b>Synonyms:</b>	RAYL; BBS19; RABL4
<b>Swissprot:</b>	Q9BW83
<b>Gene Accession:</b>	BC000566
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a GTP-binding protein that is a core component of the intraflagellar transport complex B. Characterization of the similar Chlamydomonas protein indicates a function in cell cycle control. Alternative splicing of this gene results in multiple transcript variants.