

兔抗 FNTB 多克隆抗体

中文名称：兔抗 FNTB 多克隆抗体

英文名称：Anti-FNTB rabbit polyclonal antibody

别名：farnesyltransferase, CAAX box, beta; FPTB

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FNTB

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Essential subunit of the farnesyltransferase complex . Catalyzes the transfer of a farnesyl moiety from farnesyl diphosphate to a cysteine at the fourth position from the C-terminus of several proteins having the C-terminal sequence Cys-aliphatic-aliphatic-X.
Applications:	ELISA, IHC
Name of antibody:	FNTB
Immunogen:	Fusion protein of human FNTB
Full name:	farnesyltransferase, CAAX box, beta
Synonyms:	FPTB

SwissProt:	P49356
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	40-200



