

IPO11 抗原（重组蛋白）

中文名称： IPO11 抗原（重组蛋白）

英文名称： IPO11 Antigen (Recombinant Protein)

别名： RanBP11

储存： 冷冻（-20℃）

相关类别： 抗原

概述：

Fusion protein corresponding to a region derived from 776-975 amino acids of human IPO11

技术规格：

Full name:	importin 11
Synonyms:	RanBP11
Swissprot:	Q9UI26
Gene Accession:	BC033776
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Importins, including IPO11, are a members of the karyopherin/importin-beta family of transport receptors (see KPNB1; 602738) that mediate nucleocytoplasmic transport of protein and RNA cargoes (Plafker and Macara, 2000 [PubMed 11032817]).