

LMO3 抗原（重组蛋白）

中文名称：LMO3 抗原（重组蛋白）

英文名称：LMO3 Antigen (Recombinant Protein)

别名：RBTN3; RHOM3; RBTNL2; Rhom-3

储存：冷冻（-20℃）

相关类别 抗原

概述：

Fusion protein corresponding to a region derived from 1-145 amino acids of human LMO3

技术规格：

Full name:	LIM domain only 3
Synonyms:	RBTN3; RHOM3; RBTNL2; Rhom-3
Swissprot:	Q8TAP4
Gene Accession:	BC026311
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene belongs to the rhombotin family of cysteine-rich LIM domain oncogenes. This gene is predominantly expressed in the brain. Related family members, LMO1 and LMO2 on chromosome 11, have been reported to be involved in chromosomal translocations in T-cell leukemia. Many alternatively spliced transcript variants have been found for this gene.