

LYRM1 抗原（重组蛋白）

中文名称：LYRM1 抗原（重组蛋白）

英文名称：LYRM1 Antigen (Recombinant Protein)

别名：A211C6.1

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-122 amino acids of human LYRM1

技术规格：

Full name:	LYR motif containing 1
Synonyms:	A211C6.1
Swissprot:	O43325
Gene Accession:	BC017039
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene belongs to the mitochondrial leucine/tyrosine/arginine motif family of proteins. Proteins of this family are short polypeptides that contain a leucine/tyrosine/arginine motif near the N-terminus. This gene is widely expressed with high levels in omental adipose tissue of obese individuals. In adipose tissue, the protein is localized to the nucleus where it promotes preadipocyte proliferation and lowers the rate of apoptosis to regulate adipogenesis.

se tissue homeostasis. Overexpression of this gene in adipocytes causes abnormal mitochondrial morphology and mitochondrial dysfunction. Alternative splicing results in multiple transcript variants.