

## MELTF 抗原（重组蛋白）

中文名称：MELTF 抗原（重组蛋白）

英文名称：MELTF Antigen (Recombinant Protein)

别名：melanotransferrin; MTF; MFI2; MTF1; CD228; MAP97

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 20-302 amino acids of human MELTF

技术规格：

<b>Full name:</b>	melanotransferrin
<b>Synonyms:</b>	MTf; MFI2; MTF1; CD228; MAP97
<b>Swissprot:</b>	P08582
<b>Gene Accession:</b>	BC001875
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene is a cell-surface glycoprotein found on melanoma cells. The protein shares sequence similarity and iron-binding properties with members of the transferrin superfamily. The importance of the iron binding function has not yet been identified. This gene resides in the same region of chromosome 3 as members of the transferrin superfamily. Alternative splicing results in two transcript variants.