

兔抗 FOXO1(Ab-256)多克隆抗体

- 中文名称： 兔抗 FOXO1(Ab-256)多克隆抗体
- 英文名称： Anti-FOXO1(Ab-256) rabbit polyclonal antibody
- 别名： forkhead box O1; FKH1; FKHR; FOXO1A
- 相关类别： 一抗
- 抗原： FOXO1(Ab-256)
- 储存： 冷冻（-20℃）
- 宿主： Rabbit
- 反应种属： Human, Mouse, Rat
- 标记物： Unconjugate
- 克隆类型： rabbit polyclonal

技术规格

Background:	This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in myogenic growth and differentiation. Translocation of this gene with PAX3 has been associated with alveolar rhabdomyosarcoma.
Applications:	WB, IHC, IF
Name of antibody:	FOXO1(Ab-256)
Immunogen:	Synthetic peptide of human FOXO1
Full name:	forkhead box O1

Synonyms:	FKH1; FKHR; FOXO1A
SwissProt:	Q12778
IHC positive control:	Human breast carcinoma
IHC Recommend dilution:	50-100
WB Predicted band size:	70 kDa
WB Positive control:	Hela cells
WB Recommended dilution:	500-1000
IF positive control:	Hela cells
IF Recommend dilution:	100-200



