

兔抗 FRAT1 多克隆抗体

中文名称: 兔抗 FRAT1 多克隆抗体

英文名称: Anti-FRAT1 rabbit polyclonal antibody

别名: FRAT1, WNT signaling pathway regulator

相关类别: 一抗

抗原: FRAT1

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the GSK-3-binding protein family. The protein inhibits GSK-3-mediated phosphorylation of beta-catenin and positively regulates the Wnt signaling pathway. It may function in tumor progression and in lymphomagenesis.
Applications:	ELISA, IHC
Name of antibody:	FRAT1
Immunogen:	Synthetic peptide of human FRAT1
Full name:	FRAT1, WNT signaling pathway regulator
SwissProt:	Q92837
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human tonsil
IHC Recommend dilution:	25-50

