

## 兔抗 FRMD5 多克隆抗体

中文名称：兔抗 FRMD5 多克隆抗体

英文名称：Anti-FRMD5 rabbit polyclonal antibody

别名：FERM domain containing 5

相关类别：一抗

抗原：FRMD5

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	May be involved in regulation of cell migration (PubMed: 22846708, PubMed:25448675). May regulate cell-matrix interactions via its interaction with ITGB5 and modifying ITGB5 cytoplasmic tail interactions such as with FERMT2 and TLN1. May regulate ROCK1 kinase activity possibly involved in regulation of actin stress fiber formation (PubMed:25448675).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	FRMD5
<b>Immunogen:</b>	Fusion protein of human FRMD5
<b>Full name:</b>	FERM domain containing 5

<b>SwissProt:</b>	Q7Z6J6
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	200-400

