

兔抗 FSHR 多克隆抗体

中文名称：兔抗 FSHR 多克隆抗体

英文名称： Anti-FSHR rabbit polyclonal antibody

别名： LGR1, ODG1, FSHRO

抗原： FSHR

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to family 1 of G-protein coupled receptors. It is the receptor for follicle stimulating hormone and functions in gonad development. Mutations in this gene cause ovarian dysgenesis type 1, and also ovarian hyperstimulation syndrome. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	FSHR
Immunogen:	Fusion protein of human FSHR
Full name:	follicle stimulating hormone receptor
Synonyms :	LGR1, ODG1, FSHRO

SwissProt:	P23945
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human breast cancer and Human liver cancer
IHC Recommend dilution:	50-200

