

兔抗 FST 多克隆抗体

中文名称：兔抗 FST 多克隆抗体

英文名称：Anti-FST rabbit polyclonal antibody

别名：FS

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：FST

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Binds directly to activin and functions as an activin antagonist. Specific inhibitor of the biosynthesis and secretion of pituitary follicle stimulating hormone (FSH).
Applications:	WB
Name of antibody:	FST
Immunogen:	Synthesized peptide derived from internal of human FST.
Full name:	follistatin
Synonyms :	FS
SwissProt:	P19883
WB Predicted band size:	38 kDa
WB Positive control:	HepG2, K562, COLO205, Jurkat and Hela cells lysates

WB Recommended dilution: 500-3000



