

## 兔抗 FTO 多克隆抗体

中文名称：兔抗 FTO 多克隆抗体

英文名称： Anti-FTO rabbit polyclonal antibody

别名： ALKBH9

相关类别： 一抗

抗原： FTO

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a nuclear protein of the AlkB related non-haem iron and 2-oxoglutarate-dependent oxygenase superfamily but the exact physiological function of this gene is not known. Other non-heme iron enzymes function to reverse alkylated DNA and RNA damage by oxidative demethylation. Studies in mice and humans indicate a role in nervous and cardiovascular systems and a strong association with body mass index, obesity risk, and type 2 diabetes.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	FTO

<b>Immunogen:</b>	Fusion protein of human FTO
<b>Full name:</b>	fat mass and obesity associated
<b>Synonyms:</b>	ALKBH9
<b>SwissProt:</b>	Q9C0B1
<b>IHC positive control:</b>	Human gastric cancer and human liver cancer
<b>IHC Recommend dilution:</b>	100-300
<b>WB Predicted band size:</b>	58 kDa
<b>WB Positive control:</b>	Human fetal brain tissue
<b>WB Recommended dilution:</b>	500-2000



