

兔抗 FUBP3 多克隆抗体

中文名称：兔抗 FUBP3 多克隆抗体

英文名称：Anti-FUBP3 rabbit polyclonal antibody

别名：far upstream element binding protein 3; FBP3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FUBP3

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May interact with single-stranded DNA from the far-upstream element (FUSE). May activate gene expression.
Applications:	ELISA, IHC
Name of antibody:	FUBP3
Immunogen:	Fusion protein of human FUBP3
Full name:	far upstream element binding protein 3
Synonyms:	FBP3
SwissProt:	Q96I24
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and Human tonsil
IHC Recommend dilution:	100-300



