

兔抗 GLIS1 多克隆抗体

中文名称: 兔抗 GLIS1 多克隆抗体

英文名称: Anti-GLIS1 rabbit polyclonal antibody

别名: GLIS family zinc finger 1

抗原: GLIS1

储存: 冷冻 (-20°C)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	GLIS1 is a GLI (MIM 165220)-related Kruppel-like zinc finger protein that functions as an activator and repressor of transcription (Kim et al., 2002 [PubMed 12042312]). Binds to the consensus sequence 5'-GACCACCCAC-3' (By similarity).
Applications:	ELISA, IHC
Name of antibody:	GLIS1
Immunogen:	Synthetic peptide of human GLIS1
Full name:	GLIS family zinc finger 1
SwissProt:	Q8NBF1
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human colorectal cancer and Human cervical cancer
IHC Recommend dilution:	50-300



