

兔抗 GLS2 多克隆抗体

中文名称: 兔抗 GLS2 多克隆抗体

英文名称: Anti-GLS2 rabbit polyclonal antibody

别名: glutaminase 2; GA; GLS; LGA; hLGA

抗原: GLS2

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a mitochondrial phosphate-activated glutaminase that catalyzes the hydrolysis of glutamine to stoichiometric amounts of glutamate and ammonia. Originally thought to be liver-specific, this protein has been found in other tissues as well. Alternative splicing results in multiple transcript variants that encode different isoforms.
Applications:	ELISA, WB, IHC
Name of antibody:	GLS2
Immunogen:	Synthetic peptide of human GLS2
Full name:	glutaminase 2

Synonyms:	GA; GLS; LGA; hLGA
SwissProt:	Q9UI32
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer
IHC Recommend dilution:	40-200
WB Predicted band size:	66 kDa
WB Positive control:	Mouse liver tissue lysate
WB Recommended dilution:	500-2000



