

## 兔抗 GLT8D1 多克隆抗体

中文名称：兔抗 GLT8D1 多克隆抗体

英文名称：Anti-GLT8D1 rabbit polyclonal antibody

别名：glycosyltransferase 8 domain containing 1; AD-017; MSTP139

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：GLT8D1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the glycosyltransferase family. The specific function of this protein has not been determined. Alternative splicing results in multiple transcript variants of this gene
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	GLT8D1
<b>Immunogen:</b>	Fusion protein of human GLT8D1
<b>Full name:</b>	glycosyltransferase 8 domain containing 1
<b>Synonyms:</b>	AD-017; MSTP139
<b>SwissProt:</b>	Q68CQ7
<b>ELISA Recommended dilution:</b>	2000-5000

<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	100-300
<b>WB Predicted band size:</b>	42 kDa
<b>WB Positive control:</b>	293T, LO2 cell lysates
<b>WB Recommended dilution:</b>	500-2000



