

PSMA3 抗原（重组蛋白）

中文名称：PSMA3 抗原（重组蛋白）

英文名称：PSMA3 Antigen (Recombinant Protein)

别名：proteasome subunit alpha 3; HC8; PSC3

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 2-248 amino acids of human PSMA3

技术规格：

Full name:	proteasome subunit alpha 3
Synonyms:	HC8; PSC3
Swissprot:	P25788
Gene Accession:	BC005265
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The proteasome is a multicatalytic proteinase complex with a highly ordered ring-shaped 20S core structure. The core structure is composed of 4 rings of 28 non-identical subunits; 2 rings are composed of 7 alpha subunits and 2 rings are composed of 7 beta subunits. Proteasomes are distributed throughout eukaryotic cells at a high concentration and cleave peptides in an ATP/ubiquitin-dependent process in a non-lysosomal pathway. An essential function

n of a modified proteasome, the immunoproteasome, is the processing of class I MHC peptides. This gene encodes a member of the peptidase T1A family, that is a 20S core alpha subunit. Two alternative transcripts encoding different isoforms have been identified.