

PSMB10 抗原（重组蛋白）

中文名称：PSMB10 抗原（重组蛋白）

英文名称：PSMB10 Antigen (Recombinant Protein)

别名：proteasome subunit beta 10; LMP10; MECL1; beta2i

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 40-273 amino acids of human PSMB10

技术规格：

| | |
|---------------------------|--|
| Full name: | proteasome subunit beta 10 |
| Synonyms: | LMP10; MECL1; beta2i |
| Swissprot: | P40306 |
| Gene Accession: | BC017198 |
| Purity: | >85%, as determined by Coomassie blue stained SDS-PAGE |
| Expression system: | Escherichia coli |
| Tags: | His tag C-Terminus, GST tag N-Terminus |
| Background: | The proteasome is a multicatalytic proteinase complex with a highly ordered ring-shaped 20S core structure. The core structure is composed of 4 rings of 28 non-identical subunits; 2 rings are composed of 7 alpha subunits and 2 rings are composed of 7 beta subunits. Proteasomes are distributed throughout eukaryotic cells at a high concentration and cleave peptides in an ATP/ubiquitin-dependent process in a non-lysosomal pathway. An essential function of |

a modified proteasome, the immunoproteasome, is the processing of class I MHC peptides. This gene encodes a member of the proteasome B-type family, also known as the T1B family, that is a 20S core beta subunit. Proteolytic processing is required to generate a mature subunit. Expression of this gene is induced by gamma interferon, and this gene product replaces catalytic subunit 2 (proteasome beta 7 subunit) in the immunoproteasome.