

## PTGES3L-AARSD1 抗原（重组蛋白）

中文名称：PTGES3L-AARSD1 抗原（重组蛋白）

英文名称：PTGES3L-AARSD1 Antigen (Recombinant Protein)

别名：PTGES3L-AARSD1 readthrough

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 213-412 amino acids of human PTGES3L-AARSD1

技术规格：

<b>Full name:</b>	PTGES3L-AARSD1 readthrough
<b>Swissprot:</b>	Q9BTE6
<b>Gene Accession:</b>	BC004172
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This locus represents naturally occurring readthrough transcription between the neighboring PTGES3L (prostaglandin E synthase 3 (cytosolic)-like) and AARSD1(alanyl-tRNA synthetase domain containing 1) genes on chromosome 17. The readthrough transcript encodes a fusion protein that shares sequence identity with each individual gene product.