

## RABGGTA 抗原（重组蛋白）

中文名称：RABGGTA 抗原（重组蛋白）

英文名称：RABGGTA Antigen (Recombinant Protein)

别名：Rab geranylgeranyltransferase alpha subunit; PTAR3

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 368-567 amino acids of human RABGGTA

技术规格：

<b>Full name:</b>	Rab geranylgeranyltransferase alpha subunit
<b>Synonyms:</b>	PTAR3
<b>Swissprot:</b>	Q92696
<b>Gene Accession:</b>	BC003093
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Catalyzes the transfer of a geranylgeranyl moiety from geranylgeranyl diphosphate to both cysteines of Rab proteins with the C-terminal sequence -XXCC, -XCXC and -CCXX, such as RAB1A, RAB3A, RAB5A and RAB7A.