

RBM15 抗原（重组蛋白）

中文名称：RBM15 抗原（重组蛋白）

英文名称：RBM15 Antigen (Recombinant Protein)

别名：OTT; OTT1; SPEN

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 778-977 amino acids of human RBM15

技术规格：

Full name:	RNA binding motif protein 15
Synonyms:	OTT; OTT1; SPEN
Swissprot:	Q96T37
Gene Accession:	BC006397
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Members of the SPEN (Split-end) family of proteins, including RBM15, have repressor function in several signaling pathways and may bind to RNA through interaction with spliceosome components (Hiriart et al., 2005 [PubMed 16129689]).