

RCOR3 抗原（重组蛋白）

中文名称：RCOR3 抗原（重组蛋白）

英文名称：RCOR3 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述：

Fusion protein corresponding to a region derived from 59-258 amino acids of human RCOR3

技术规格：

Full name:	REST corepressor 3
Swissprot:	Q9P2K3
Gene Accession:	BC031608
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May act as a component of a corepressor complex that represses transcription.