

RIMS3 抗原（重组蛋白）

中文名称：RIMS3 抗原（重组蛋白）

英文名称： RIMS3 Antigen (Recombinant Protein)

别名： regulating synaptic membrane exocytosis 3; NIM3; RIM3; RIM 3

储存： 冷冻（-20℃）

相关类别： 抗原

概述：

Fusion protein corresponding to a region derived from 1-308 amino acids of human RIMS3

技术规格：

Full name:	regulating synaptic membrane exocytosis 3
Synonyms:	NIM3; RIM3; RIM 3
Swissprot:	Q9UJD0
Gene Accession:	BC003103
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Regulates synaptic membrane exocytosis.