

RNF148 抗原（重组蛋白）

中文名称：RNF148 抗原（重组蛋白）

英文名称：RNF148 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述：

Fusion protein corresponding to a region derived from 35-305 amino acids of human RNF148

技术规格：

Full name:	ring finger protein 148
Swissprot:	Q8N7C7
Gene Accession:	BC029264
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	RNF148 is a single-pass membrane protein. It contains 1 PA (prot ease associated) domain and 1 RING-type zinc finger. The exact f unction of RNF148 remains unknown.