

SAMSN1 抗原（重组蛋白）

中文名称：SAMSN1 抗原（重组蛋白）

英文名称：SAMSN1 Antigen (Recombinant Protein)

别名：SLy2; HACs1; NASH1; SASH2; SH3D6B

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 174-373 amino acids of human SAMSN1

技术规格：

Full name:	SAM domain, SH3 domain and nuclear localization signals 1
Synonyms:	SLy2; HACs1; NASH1; SASH2; SH3D6B
Swissprot:	Q9NSI8
Gene Accession:	BC029112
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	SAMSN1 is a member of a novel gene family of putative adaptors and scaffold proteins containing SH3 and SAM (sterile alpha motif) domains (Claudio et al., 2001 [PubMed 11536050]).