

## SCAND1 抗原（重组蛋白）

中文名称：SCAND1 抗原（重组蛋白）

英文名称： SCAND1 Antigen (Recombinant Protein)

别名： RAZ1; SDP1

储存： 冷冻（-20℃）

相关类别 抗原

概述：

Fusion protein corresponding to a region derived from 1-179 amino acids of human SCAND1

技术规格：

<b>Full name:</b>	SCAN domain containing 1
<b>Synonyms:</b>	RAZ1; SDP1
<b>Swissprot:</b>	P57086
<b>Gene Accession:</b>	BC000785
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a SCAN box domain-containing protein. The SCAN domain is a highly conserved, leucine-rich motif of approximately 60 aa originally found within a subfamily of zinc finger proteins. This gene belongs to a family of genes that encode an isolated SCAN domain, but no zinc finger motif. This protein binds to and may regulate the function of the transcription factor myeloid zinc finger 1B. Alternate splicing results i

n multiple transcript variants.