

SCGN 抗原（重组蛋白）

中文名称：SCGN 抗原（重组蛋白）

英文名称：SCGN Antigen (Recombinant Protein)

别名：secretagoin, EF-hand calcium binding protein; SEGN; CALBL; SECRET; setagin; DJ501N12.8

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-276 amino acids of human SCGN

技术规格：

Full name:	secretagoin, EF-hand calcium binding protein
Synonyms:	SEGN; CALBL; SECRET; setagin; DJ501N12.8
Swissprot:	O76038
Gene Accession:	BC000336
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The encoded protein is a secreted calcium-binding protein which is found in the cytoplasm. It is related to calbindin D-28K and c alretinin. This protein is thought to be involved in KCL-stimulated calcium flux and cell proliferation.