

## SCNM1 抗原（重组蛋白）

中文名称：SCNM1 抗原（重组蛋白）

英文名称：SCNM1 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述：

Fusion protein corresponding to a region derived from 1-230 amino acids of human SCNM1

技术规格：

<b>Full name:</b>	sodium channel modifier 1
<b>Swissprot:</b>	Q9BWG6
<b>Gene Accession:</b>	BC000264
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	SCNM1 is a zinc finger protein and putative splicing factor. In mice, Scnm1 modifies phenotypic expression of Scn8a (MIM 600702) mutations (Buchner et al., 2003 [PubMed 12920299]). Plays a role in alternative splicing of pre-mRNAs, possibly by contributing to the selection of non-consensus donor sites.