

## SCRN2 抗原（重组蛋白）

中文名称：SCRN2 抗原（重组蛋白）

英文名称：SCRN2 Antigen (Recombinant Protein)

别名：secernin 2; Ses2

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-200 amino acids of human SCRN2

技术规格：

<b>Full name:</b>	secernin 2
<b>Synonyms:</b>	Ses2
<b>Swissprot:</b>	Q96FV2
<b>Gene Accession:</b>	BC017317
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	SCRN genes encode protein members of the secernin family including Secernin 1, Secernin 2 and Secernin 3. Function of Secernin 2 and Secernin 3 are so far not well understood. Secernin 1 is a novel 50-kDa cytosolic protein that appears to be involved in the regulation of exocytosis from peritoneal mast cells. Secernin-1 a novel tumor-associated antigen (TAA) and may be a universal marker of different cancer types including gastric cancer, which has been validated by several

researches. Measurement of SCRN1 can be used in the early detection of cancer or in the surveillance of patients who undergo surgery.