

兔抗 GNAI1 多克隆抗体

中文名称：兔抗 GNAI1 多克隆抗体

英文名称： Anti-GNAI1 rabbit polyclonal antibody

别名： G protein subunit alpha i1; Gi

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： GNAI1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Guanine nucleotide binding proteins are heterotrimeric signal-transducing molecules consisting of alpha, beta, and gamma subunits. The alpha subunit binds guanine nucleotide, can hydrolyze GTP, and can interact with other proteins. The protein encoded by this gene represents the alpha subunit of an inhibitory complex. The encoded protein is part of a complex that responds to beta-adrenergic signals by inhibiting adenylate cyclase. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB, IHC

Name of antibody:	GNAI1
Immunogen:	Fusion protein of human GNAI1
Full name:	G protein subunit alpha i1
Synonyms:	Gi
SwissProt:	P63096
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer
IHC Recommend dilution:	100-200
WB Predicted band size:	40 kDa
WB Positive control:	A172 and Hela cell lysates
WB Recommended dilution:	1000-5000



