

SHC4 抗原（重组蛋白）

中文名称：SHC4 抗原（重组蛋白）

英文名称：SHC4 Antigen (Recombinant Protein)

别名：RaLP; SHCD

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 431-630 amino acids of human SHC4

技术规格：

Full name:	SHC adaptor protein 4
Synonyms:	RaLP; SHCD
Swissprot:	Q6S5L8
Gene Accession:	BC033907
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Activates both Ras-dependent and Ras-independent migratory pathways in melanomas. Contributes to the early phases of agrin-induced tyrosine phosphorylation of CHRN1.