

## SIGLEC6 抗原（重组蛋白）

中文名称：SIGLEC6 抗原（重组蛋白）

英文名称：SIGLEC6 Antigen (Recombinant Protein)

别名：CD327; CD33L; OBBP1; CD33L1; CD33L2; CDW327

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 254-453 amino acids of human SIGLEC6

技术规格：

<b>Full name:</b>	sialic acid binding Ig like lectin 6
<b>Synonyms:</b>	CD327; CD33L; OBBP1; CD33L1; CD33L2; CDW327
<b>Swissprot:</b>	O43699
<b>Gene Accession:</b>	BC035359
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a member of the SIGLEC (sialic acid binding immunoglobulin-like lectin) family of proteins. The encoded transmembrane receptor binds sialyl-TN glycans and leptin. Placental expression of the encoded protein is upregulated in preeclampsia.