

## 兔抗 GNG5 多克隆抗体

中文名称：兔抗 GNG5 多克隆抗体

英文名称： Anti-GNG5 rabbit polyclonal antibody

抗 原： GNG5

储 存： 冷冻（-20℃） 避光

宿 主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	G proteins are trimeric (alpha-beta-gamma) membrane-associated proteins that regulate flow of information from cell surface receptors to a variety of internal metabolic effectors. Interaction of a G protein with its activated receptor promotes exchange of GTP for GDP that is bound to the alpha subunit.
<b>Applications:</b>	IF
<b>Name of antibody:</b>	GNG5
<b>Immunogen:</b>	Synthesized peptide derived from internal of human GNG5.
<b>Full name:</b>	guanine nucleotide binding protein (G protein), gamma 5
<b>Synonyms :</b>	0

<b>SwissProt:</b>	P63218
<b>IF positive control:</b>	A549 cells
<b>IF Recommend dilution:</b>	100-500

