

SMOC1 抗原（重组蛋白）

中文名称：SMOC1 抗原（重组蛋白）

英文名称：SMOC1 Antigen (Recombinant Protein)

别名：SPARC related modular calcium binding 1; OAS

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 236-435 amino acids of human SMOC1

技术规格：

Full name:	SPARC related modular calcium binding 1
Synonyms:	OAS
Swissprot:	Q9H4F8
Gene Accession:	BC011548
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a multi-domain secreted protein that may have a critical role in ocular and limb development. Mutations in this gene are associated with microphthalmia and limb anomalies. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.