

兔抗 GNL3L 多克隆抗体

中文名称：兔抗 GNL3L 多克隆抗体

英文名称： Anti-GNL3L rabbit polyclonal antibody

抗 原： GNL3L

储 存： 冷冻（-20℃） 避光

宿 主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

| | |
|--------------------------------|--|
| Background: | The protein encoded by this gene appears to be a nucleolar GTPase that is essential for ribosomal pre-rRNA processing and cell proliferation. Two transcript variants encoding the same protein have been found for this gene. |
| Applications: | WB |
| Name of antibody: | GNL3L |
| Immunogen: | Synthesized peptide derived from internal of human GNL3L. |
| Full name: | guanine nucleotide binding protein-like 3 (nucleolar)-like |
| Synonyms : | 0 |
| SwissProt: | Q9NVN8 |
| WB Predicted band size: | 66 kDa |
| WB Positive control: | Hela cells, 293 cells and HT-29 cells lysates |

WB Recommended dilution: 500-3000

