

SNAPC1 抗原（重组蛋白）

中文名称：SNAPC1 抗原（重组蛋白）

英文名称：SNAPC1 Antigen (Recombinant Protein)

别名：small nuclear RNA activating complex polypeptide 1; SNAP43; PTFgamma

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 169-368 amino acids of human SNAPC1

技术规格：

Full name:	small nuclear RNA activating complex polypeptide 1
Synonyms:	SNAP43; PTFgamma
Swissprot:	Q16533
Gene Accession:	BC014984
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Part of the SNAPc complex required for the transcription of both RNA polymerase II and III small-nuclear RNA genes. Binds to the proximal sequence element (PSE), a non-TATA-box basal promoter element common to these 2 types of genes. Recruits TBP and BRF2 to the U6 snRNA TATA box.