

## 兔抗 GNPAT 多克隆抗体

中文名称: 兔抗 GNPAT 多克隆抗体

英文名称: Anti-GNPAT rabbit polyclonal antibody

别名: DAPAT; DAP-AT; DHAPAT

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: GNPAT

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes an enzyme located in the peroxisomal membrane which is essential to the synthesis of ether phospholipids. Mutations in this gene are associated with rhizomelic chondrodysplasia punctata.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	GNPAT
<b>Immunogen:</b>	Synthetic peptide of human GNPAT
<b>Full name:</b>	glyceronephosphate O-acyltransferase
<b>Synonyms :</b>	DAPAT; DAP-AT; DHAPAT
<b>SwissProt:</b>	O15228

<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human cervical cancer
<b>IHC Recommend dilution:</b>	50-200

