

## SNRPB2 抗原（重组蛋白）

中文名称：SNRPB2 抗原（重组蛋白）

英文名称：SNRPB2 Antigen (Recombinant Protein)

别名：Msl1; U2B''

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 26-225 amino acids of human SNRPB2

技术规格：

<b>Full name:</b>	small nuclear ribonucleoprotein polypeptide B2
<b>Synonyms:</b>	Msl1; U2B''
<b>Swissprot:</b>	P08579
<b>Gene Accession:</b>	BC018022
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene associates with stem loop IV of U2 small nuclear ribonucleoprotein (U2 snRNP) in the presence of snRNP-A'. The encoded protein may play a role in pre-mRNA splicing. Autoantibodies from patients with systemic lupus erythematosus frequently recognize epitopes on the encoded protein. Two transcript variants encoding the same protein have been found for this gene.