

STK32C 抗原（重组蛋白）

中文名称：STK32C 抗原（重组蛋白）

英文名称：STK32C Antigen (Recombinant Protein)

别名：serine/threonine kinase 32C; PKE; YANK3

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 287-486 amino acids of human STK32C

技术规格：

Full name:	serine/threonine kinase 32C
Synonyms:	PKE; YANK3
Swissprot:	Q86UX6
Gene Accession:	BC015792
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a member of the serine/threonine protein kinase family. It is thought that this family member is functional in brain due to its high expression levels there. DNA methylation differences have been found in this gene in monozygotic twins that are discordant for adolescent depression. Alternative splicing results in multiple transcript variants.