

兔抗 GOLGA7 多克隆抗体

中文名称: 兔抗 GOLGA7 多克隆抗体

英文名称: Anti-GOLGA7 rabbit polyclonal antibody

别名: golgin A7; GCP16; GOLGA7A; HSPC041; GOLGA3AP1

抗原: GOLGA7

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	May be involved in protein transport from Golgi to cell surface. The ZDHHC9-GOLGA7 complex is a palmitoyltransferase specific for HRAS and NRAS.
Applications:	ELISA, IHC
Name of antibody:	GOLGA7
Immunogen:	Fusion protein of human GOLGA7
Full name:	golgin A7
Synonyms:	GCP16; GOLGA7A; HSPC041; GOLGA3AP1
SwissProt:	Q7Z5G4
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer

IHC Recommend dilution: 40-200



