

SULT4A1 抗原（重组蛋白）

中文名称：SULT4A1 抗原（重组蛋白）

英文名称：SULT4A1 Antigen (Recombinant Protein)

别名：NST; BRSTL1; SULTX3; BR-STL-1; DJ388M5.3; hBR-STL-1

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-260 amino acids of human SULT4A1

技术规格：

Full name:	sulfotransferase family 4A member 1
Synonyms:	NST; BRSTL1; SULTX3; BR-STL-1; DJ388M5.3; hBR-STL-1
Swissprot:	Q9BR01
Gene Accession:	BC028171
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the sulfotransferase family. The encoded protein is a brain-specific sulfotransferase believed to be involved in the metabolism of neurotransmitters. Polymorphisms in this gene may be associated with susceptibility to schizophrenia.