

兔抗 GOLIM4 多克隆抗体

中文名称: 兔抗 GOLIM4 多克隆抗体

英文名称: Anti-GOLIM4 rabbit polyclonal antibody

别名: golgi integral membrane protein 4; P138; GIMPC; GOLPH4; GPP130

抗原: GOLIM4

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The Golgi complex plays a key role in the sorting and modification of proteins exported from the endoplasmic reticulum. The protein encoded by this gene is a type II Golgi-resident protein. It may process proteins synthesized in the rough endoplasmic reticulum and assist in the transport of protein cargo through the Golgi apparatus.
Applications:	ELISA, IHC
Name of antibody:	GOLIM4
Immunogen:	Synthetic peptide of human GOLIM4
Full name:	golgi integral membrane protein 4

Synonyms:	P138; GIMPC; GOLPH4; GPP130
SwissProt:	O00461
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer
IHC Recommend dilution:	20-100

