

## TEAD3 抗原（重组蛋白）

中文名称：TEAD3 抗原（重组蛋白）

英文名称：TEAD3 Antigen (Recombinant Protein)

别名：TEF5; TEAD5; TEF-5; DTEF-1; ETFR-1; TEAD-3

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 236-435 amino acids of human TEAD3

技术规格：

<b>Full name:</b>	TEA domain transcription factor 3
<b>Synonyms:</b>	TEF5; TEAD5; TEF-5; DTEF-1; ETFR-1; TEAD-3
<b>Swissprot:</b>	Q99594
<b>Gene Accession:</b>	BC027877
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene product is a member of the transcriptional enhancer factor (TEF) family of transcription factors, which contain the TEA/ATTS DNA-binding domain. It is predominantly expressed in the placenta and is involved in the transactivation of the chorionic somatomammotropin-B gene enhancer. Translation of this protein is initiated at a non-AUG (AUA) start codon.